

# SURGICAL REGION

Manufacture & Exporters




## Eye Instruments


 [www.surgicalregion.com](http://www.surgicalregion.com)



<b>Description</b>	<b>Page</b>
EYE SPECULA	2-8
SCALPELS AND KNIVES	10,11
KNIVES	13-15
HOOKS	16
RETRACTORS	17
LOOPS	18
SPATULAS	19-21
CHALAZION CURETTES	22,24
FOREIGN BODY INSTRUMENTS	23
LENS MANUPILATORS	25,26
IRIS HOOKS	27
CAPSULE INSTRUMENTS	28
DISSECTING AND TISSUE FORCEPS	30,31
IRIS FORCEPS	32-34
FIXATION FORCEPS	35-39
TYING FORCEPS	40-44
SUTURE FORCEPS	45-49
CORNEAL FORCEPS	50-56
CAPSULE FORCEPS	57
CAPSULOREHEXIS FORCEPS	58
NUCLEUS FORCEPS	59
FRAGMENT FORCEPS	60
CHALAZION FORCEPS	61,62
MUSCLE FORCEPS	63
MICRO NEEDLE HOLDERS	65-69
DELICATE SCISSORS	71-75
MICRO SCISSORS	76-87
LACRIMAL DILATOR	89,90
DCR-INSTRUMENTS	91
CANNULA	92-96

- Eye surgery
- Augen-Instrumentarium
- Oftalmologia
- Chirurgie ophtalmologique
- Oftalmologia






**Schott Eye Speculum**  
Fenestrated, individually adjustable blades for controlled lid elevation  
**TAB-100.01** Adult Size  
**TAB-100.02** Child Size




**Schott Eye Speculum**  
With solid, individually adjustable blades for controlled lid elevation  
**TAB-102.01** Adult Size  
**TAB-102.02** Child Size


**Weiss Eye Speculum**  
With open blades  
**PAB-104.06**




**Weiss Eye Specula Large**  
With fenestrated blades  
**PAB-105.07**



**Weiss Eye Specula small**  
With fenestrated blades  
**PAB-106.06**




**Weiss Eye Specula large**  
With solid blades  
**PAB-107.09**




**Weiss Eye Specula small**  
With solid blades  
**PAB-108.08**



**Lester-Burch Eye Specula**  
With locking screw, solid blades  
PAB-109.14




**Lancaster Eye Specula**  
With locking screw, solid blades  
PAB-110.06



**Castroviejo Eye Specula**  
PAB-111.14 Small (14 mm blades)  
PAB-112.16 Medium (16 mm blades)  
PAB-113.18 Large (18 mm blades)



**Mellinger-Beard Eye Specula**  
With solid blades  
PAB-114.01 Child size  
PAB-114.02 Adult size




**Mellinger Eye Speculum**

With fenestrated blades

**PAB-115.01** Child size

**PAB-116.02** Adult size



**Cook Eye Speculum**

Infant size

**PAB-117.01**



**Cook Eye Speculum**

Child size


**PAB-118.02**



**Cook Eye Speculum**

Adolescent size


**PAB-119.03**




**Cook Eye Speculum**

Adult size

**PAB-120.04**





**Barraquer Eye Speculum**


Fenestrated round blades

**PAB-121.14** Adult size (14 mm wide)

**PAB-122.10** Child size (10 mm wide)

**PAB-123.05** Infant size (5 mm wide)

**PAB-124.03** premature (3 mm wide)



**Barraquer Eye Speculum**


Fenestrated blades

**PAB-122.15** Adult size (15 mm wide)

**PAB-122.10** Child size (10 mm wide)

**PAB-122.05** Infant size (5 mm wide)

**PAB-122.03** Premature (3 mm wide)




**Barraquer Eye Speculum**

Solid blades

**PAB-123.15** Adult size (15 mm wide)

**PAB-124.10** Child size (10 mm wide)




**Kratz - Barraquer Eye Speculum**

Open blades

**PAB-125.01** Adult size

**PAB-126.02** Child size

**PAB-127.03** Heavy 1 mm model



**Jaffe Lid Retractor**

With fenestrated blade, may be clamped to surgical drape with suture or elastic band

**PAB-128.14** Adult size (14 mm wide)

**PAB-129.10** Child size (10 mm wide)



**Tennant Lid Retractor**

With solid blade, may be clamped to surgical drape with suture or elastic band

**PAB-130.14** Adult size (14 mm wide)



**Clayman Lid Retractor (set of two)**

Mobile, fenestrated blade may be clamped to surgical drape with suture or elastic band

**PAB-131.14** Adult size (14 mm wide)



**Sauer Lid Speculum**

Self-retaining, with solid blades


**PAB-132.10** Child size (10 mm wide)

**PAB-133.06** Infant size (6 mm wide)

**PAB-134.04** Premature (4 mm wide)


**Murdoch Eye Speculum**

With large fenestrated blades  
PAB-135.01 Right model  
PAB-136.02 Left model




**Clark Eye Speculum**

With adjustment screw and open blades  
PAB-137.08




**Maumee - Park Eye Speculum**


Mobile fenestrated blades with locking screw and hooked enter blade  
PAB-138.01



**Park - Guyton Eye Speculum**

Mobile fenestrated blades with vertical suture posts and hooked center blade  
PAB-139.02






**Castroviejo Lid clamp**  
With locking screw, solid blades  
**PAB-140.09**  
**Jaeger Lid Plate**  
Stainless steel  
**PAB-141.08**



**Klar Lid Everter**  
Open  
**PAB-142.12**



**Desmarres Lid Retractor**  
With solid blade (standard pattern)  
**PAB-143.08** 8 mm wide (Conway)  
**PAB-144.10** 10 mm wide # 0  
**PAB-145.12** 12 mm wide # 1  
**PAB-146.14** 14 mm wide # 2  
**PAB-147.16** 16 mm wide # 3  
**PAB-148.18** 18 mm wide # 4



**Desmarres Lid Retractor**  
With fenestrated blade (new pattern)  
**PAB-149.12** 12 mm wide # 1  
**PAB-150.14** 14 mm wide # 2



- Eye surgery
- Augen-Instrumentarium
- Oftalmologia
- Chirurgie ophtalmologique
- Oftalmologia





**PAB-151.03**  
Scalpel Handle No. 3



**PAB-152.07**  
Scalpel Handle No. 7K



Scalpel Blade Fig. 11  
sterile, box of 100 pcs



Scalpel Blade Fig. 15  
sterile, box of 100 pcs



Disposable Scalpel No. 11 with Plastic Handle  
sterile, package of 10 pieces



Disposable Scalpel No. 15 with Plastic Handle  
sterile, package of 10 pieces



**Castroviejo Razor Blade Holder**  
Overall length 130 mm with 8 mm handle  
**PAB-153.00**



**Castroviejo Razor Blade Holder**  
With 6 mm handle  
**PAB-154.90** Overall length 90 mm  
**PAB-155.12** Overall length 120 mm



**Troutman Razor Blade Holder Straight**  
Overall length 95mm  
**PAB-156.95**




**Troutman Razor Blade Holder angled**  
Overall length 90mm  
**PAB-157.90**



**Barraquer Razor Blade Holder**  
Overall length 110 mm  
**PAB-158.11**

- Eye surgery
- Augen-Instrumentarium
- Oftalmologia
- Chirurgie ophtalmologique
- Oftalmologia





**Knapp Discission Knife**  
Overall Length 120mm  
**PAB-159.50** 5,0 mm blade  
**PAB-160.60** 6,0 mm blade  
**PAB-161.70** 7,0 mm blade



**Ziegler Discission Knife**  
Overall Length 115mm  
**PAB-162.20** 2,0 mm blade  
**PAB-163.30** 3,0 mm blade  
**PAB-164.50** 5,0 mm blade  
**PAB-165.60** 6,0 mm blade



**Wecker Discission Knife**  
Overall Length 125mm  
**PAB-166.12** Double adge



**Swiss Pattern Discission Knife**  
Overall Length 120mm  
**PAB-167.12**



**Wheeler Discission Knife**  
Overall Length 115mm  
**PAB-168.18** 18 mm cutting edge



**Barkan Goniotomy Knife**  
Overall Length 120mm  
**PAB-169.12**



**Schele Goniotomy Knife**  
Overall Length 125mm  
**PAB-170.12**



**Graefe Straight Cataract Knife**

Overall length 130mm

**PAB-171.15** 1,5 x 25 mm

**PAB-172.20** 2,0 x 30 mm

**PAB-173.23** 2,3 x 33 mm



**Graefe Cystotome Knife**

Overall length 125mm, malleable shaft with 1,5 mm cutting edge

**PAB-174.12**



**Hamburg Air Keratome**

Overall length 125mm, 12 mm, long 45dg angled blade with sharp tip and cutting edges

**PAB-175.12**



**Bonn Scleral Resection Knife**

Overall length 125mm

**PAB-176.04** 4 mm blade

**PAB-177.07** 7 mm blade



**Tooke Corneal Knife**

Overall length 115 mm, smooth blade with curved cutting edges

**PAB-178.11** 3 x 18 mm blade



**Paufique Graft Knife**

Overall length 105 mm, angled tip curved cutting edge

**PAB-179.10**



**Cataract Knife**  
 Overall length 135mm  
**PAB-180.23** Blade 32 x 2,3 mm  
**PAB-181.26** Blade 35 x 2,6 mm  
**PAB-182.30** Blade 35 x 3,0 mm



**Deuschmann Cataract Knife**  
 Overall Length 135mm  
**PAB-183.13**



**Jaeger Keratome Angled**  
 Overall length 135mm  
**PAB-184.10** 10 x 10 mm  
**PAB-185.11** 11 x 11 mm  
**PAB-186.12** 12 x 12 mm  
**PAB-187.13** 13 x 13 mm



**Jaeger Keratome Straight**  
 Overall Length 135mm  
**PAB-188.10** 10 x 10 mm  
**PAB-189.11** 11 x 11 mm  
**PAB-190.12** 12 x 12 mm  
**PAB-191.13** 13 x 13 mm



**Castroviejo Keratome Angled**  
 Overall Length 135mm  
**PAB-192.13** 4 x 13 mm



**Jameson Muscle Hook**

Overall length 140 mm

**PAB-193.01** Delicate

**PAB-194.02** Small

**PAB-195.03** Large



**Graefe Strabismus Hook**

**PAB-196.01** # -1 small

**PAB-197.02** # 2 medium

**PAB-198.03** # 3 large



**Stevens Tentotomy Hook**

Overall Length 120 mm, blunt 6mm curved tip

**PAB-199.12**



**Gass Retinal Detachment Hook**

Overall length 140 mm, flat tip with oval hole

**PAB-200.14**



**Burch Fixation Pick**

Overall length 130 mm, with ball stop

**PAB-201.13**



**Twist Fixation Hook**

Overall length 120 mm, for scleral fixation

**PAB-202.12R** Right

**PAB-203.12L** Left



**Tyrell Iris Hook**

Overall length 125 mm

**PAB-204.12** Sharp

**PAB-205.12** Blunt



**Knapp Lacrimal Sac Retractor**  
Overall length 135mm, 4 blunt prongs  
**PAB-206.13**



**Rollet Lacrimal Sac Retractor**  
Overall length 135mm, 4 sharp prongs  
**PAB-207.13**




**Arruga Globe Retractor**  
Overall length 145mm  
**PAB-208.14**



**Scherpens Orbital Retractor**  
Overall length 145mm, smooth forked blade  
**PAB-209.14**



**Snellen Lens Loop**  
Overall length 140mm  
**PAB-210.14** 15 x 4 mm smooth loop



**Wilder Lens Loop**  
Overall length 140mm  
**PAB-211.14** 15 x 4 mm serrated loop



**Keener - Arlt Lens Nucleus Loop**  
Strong curved 3mm wide serrated loop for small incision  
nucleus fragment removal  
**PAB-212.03**



**Kansas Lens Nucleus Loop**  
Light curved 3 mm wide serrated loop for small incision  
nucleus fragment removal  
**PAB-213.03**



**Lewis Lens Loop**  
**PAB-214.01** Small  
**PAB-215.02** Large




**Morrison Lens Loop and Probe**  
Very thin loop and delicate olive tip  
**PAB-216.02**




**Kirby Lens Loop and Expressor Hook**  
Overall length 145 mm  
**PAB-217.14**




**Simcoe Lens Nucleus Loop**  
Overall length 110 mm, finely serrated on top loop  
**PAB-218.11**




**Knolle Lens Nucleus Spatula**  
Overall Length 120 mm, flat 0,5 mm wide blade  
**PAB-219.12**




**Sinsky Lens Nucleus Spatula**  
Overall length 120 mm, slightly curved 30 mm long narrow blade  
**PAB-220.12**




**Girard Synechia Spatula**  
Overall length 115 mm, 0,5 mm wide 10mm vaulted blade  
**PAB-221.11**




**Jaffe Lens Spatula**  
Overall length 115mm, smooth 0,5mm wide 10 mm angled blade  
**PAB-222.11**




**Johnson Sideport Nucleus Spatula**  
Overall length 115mm, 0,25 mm wide blade, 10 mm vaulted shaft with sharp tip  
**PAB-223.11**



**Birnbaum Lens Spatula**  
Overall length 120 mm, 0,25 mm wide 10 mm vaulted blade  
**PAB-224.12**



**Hirschmann Lens Spatula**  
Overall length 115 mm, flat 0,5 mm wide 10 mm angled blade  
**PAB-225.11**



**Lippman Sideport Nucleus Spatula**  
Overall length 115 mm, 0,5 mm wide blade 10 mm vaulted shaft with blunt forked tip  
**PAB-226.11**



**Barraquer Iris Spatula**  
Overall length 120 mm, 0,5 mm round blade with smooth pointed tip  
**PAB-227.07** 7 mm angled  
**PAB-228.10** 10 mm angled  
**PAB-229.15** 15 mm angled



**Eischnig Cyclodialysis Spatula**  
Overall length 135 mm, 1,0 mm wide flat blade  
**PAB-230.13** 12 mm angled



**Castroviejo Cyclodialysis Spatula**  
Overall length 135 mm, with 1,0 mm wide flat blade  
**PAB-231.10** 10 mm angled  
**PAB-232.15** 15 mm angled



**Castroviejo Cyclodialysis Spatula**  
Overall length 135 mm, double - ended 10 and 15 mm long blade with smoothly pointed tips  
**PAB-233.50** 0,50 mm wide blade  
**PAB-234.75** 0,75 mm wide blade  
**PAB-235.99** 1,00 mm wide blade



**Paton Spatula and Spoon**  
Overall length 140mm  
**PAB-236.14**



**Wheeler Cyclodialysis Spatula**

Overall length 135mm, 1,75 mm wide blades

**PAB-237.13** Stainless steel

**PAB-238.13** Sterling silver



**Wecker Iris Spatula**

Overall length 135mm, delicate 0,5 mm wide blade

**PAB-239.13**



**Wecker Iris Spatula**

Overall length 135mm, small 1,0 mm wide blade

**PAB-240.10** Stainless steel

**PAB-241.10** Sterling silver



**Wecker Iris Spatula**

Overall length 135mm, 2,0 mm wide blade

**PAB-242.13** Stainless steel

**PAB-243.13** Sterling silver



**Kimura Platinum Spatula**

Overall length 145mm, with wide ultra thin, flat blade

**PAB-244.14**



**Sautter Orbital Spatula**

Overall length 145mm, angled concave curved blade

**PAB-245.12** 12 mm wide blade

**PAB-246.15** 15 mm wide blade

**PAB-247.18** 18 mm wide blade



**Meyhöfer Chalazion Curette**

Overall length 130 mm, sharp round cup

- |                             |                             |
|-----------------------------|-----------------------------|
| <b>PAB-248.05</b> ø 0.50 mm | <b>PAB-252.20</b> ø 2.00 mm |
| <b>PAB-249.10</b> ø 1.00 mm | <b>PAB-253.25</b> ø 2.50 mm |
| <b>PAB-250.15</b> ø 1.50 mm | <b>PAB-254.30</b> ø 3.00 mm |
| <b>PAB-251.17</b> ø 1.75 mm | <b>PAB-255.35</b> ø 3.50 mm |



**Hebra Chalazion Curette**

Overall length 130 mm, sharp oval cup

- |                             |                             |
|-----------------------------|-----------------------------|
| <b>PAB-255.10</b> ø 1.00 mm | <b>PAB-259.30</b> ø 3.00 mm |
| <b>PAB-257.20</b> ø 2.00 mm | <b>PAB-260.35</b> ø 3.50 mm |
| <b>PAB-258.25</b> ø 2.50 mm |                             |



**Axenfeld Chalazion Curette**

Overall length 145 mm double ended

**PAB-261.14**



**Bunge Enucleation Curette**

Overall length 145 mm, with round sharp spoon

- |                          |                           |
|--------------------------|---------------------------|
| <b>PAB-262.05</b> ø 6 mm | <b>PAB-264.10</b> ø 10 mm |
| <b>PAB-263.08</b> ø 8 mm | <b>PAB-265.12</b> ø 12 mm |



**Gillis- Welsh Capsule Polishing Curette**

Overall length 125 mm, 8 mm backward angled

**PAB-266.12** ø 1.5 mm



**Daviel Lens Spoon**

Overall length 135 mm, light curved, small spoon


**PAB-267.13** stainless steel



**Daviel Lens Spoon**

Overall length 135 mm, light curved, large spoon


**PAB-268.13** stainless steel




**Wells Enucleation Spoon**

Overall length 145 mm, 21 mm wide blade with notch

**PAB-269.14**




**Needle and Gouge**  
Foreign Body Spud in screw handle  
**PAB-270.00**




**Needle and Gouge**  
Foreign Body Spud in sliding handle  
**PAB-271.09**



**Ellis Foreign Body Spud**  
Overall length 120mm, very thin lightly curved  
round tip  
**PAB-272.12**




**Francis Foreign Body Spud**  
Overall length 120mm, flag shaped angled  
tip with sharp edge  
**PAB-273.12**





**Dezel Foreign Body Needle**  
Overall length 120mm, sharp spear tip  
**PAB-274.12** Straight  
**PAB-275.12** Curved




**Golf Club Foreign Body Spud**  
Overall length 115mm, angled lightly rounded  
with sharp edge  
**PAB-276-11**




**Kansas Nucleus Vectis**  
 Overall length 115 mm, cutting board for  
 Kansas nucleus bisector and trisector  
**PAB-280.01** For right handed surgeons  
**PAB-280.02** For left handed surgeons



**Kansas Nucleus Bisector**  
 Overall length 110 mm, 10 mm angled blade with semi - sharp cutting edge  
 on the vertical plane, to divide the nucleus in two parts  
**PAB-281.01**




**Kansas Nucleus Trisector**  
 Overall length 110 mm  
**PAB-282.01** For right handed surgeons  
**PAB-282.02** For left handed surgeons




**Clayman Nucleus Rotator**  
 Overall length 115 mm, V- Notched 10 mm angled tip  
**PAB-284.11** Flat 0,5 mm Wide blade





**Jaffe-Bechert Lens Nucleus Rotator**  
 Overall length 115 mm, 9 mm angled  
**PAB-286.11** Blunt forked tip




**Osher Nucleus Manipulator**  
 Overall length 125 mm, with 0,5 mm two - prong tip  
**PAB-287.12** Straight  
**PAB-287.22** Angled




**Kansas Nucleus Loop**  
 Overall length 135 mm, light curved 3 mm wide serrated loop  
**PAB-288.13**




**Gruber - Hofmann Irrigating Vectis**  
 25 ga x 3 mm wide serrated loop with two front ports  
**PAB-290.25**




**Fenzl Lens Manipulating Hook**  
Overall length 120 mm 0,15 mm blunt T-Shaped tip,  
**PAB-292.01** Straight  
**PAB-292.02** Angled




**Lester Lens Pusher**  
Overall length 115 mm straight shaft with Ø 0,25 mm,  
hour-glass shaped tip  
**PAB-293.11**




**Lester Lens Pusher**  
Overall length 115 mm angled shaft with Ø 0,25 mm  
hour-glass shaped tip  
**PAB-294.11**




**Osher Lens Manipulator**  
Overall length 115 mm, Y-shaped Hook  
**PAB-296.11** Straight  
**PAB-297.11** Angled




**Jaffe IOL Manipulating Hook**  
Overall length 125mm, blunt tip with 0,15 mm  
**PAB-298.12** Straight  
**PAB-299.12** Angled



**Sinskey II IOL Manipulating Hook**  
Overall length 125mm, blunt tip 0,2 mm  
**PAB-300.12** Straight  
**PAB-301.12** Angled




**Sinsky II IOL Manipulating Hook**  
Overall length 115 mm, blunt tip  
**PAB-302.11** Straight 0,2 mm




**Lewicky IOL Manipulating Hook**  
Overall length 125 mm, blunt tip  
**PAB-304.12** 0,15 mm




**Rentsch Manipulating Hook and Iris Retractor**  
Overall length 115 mm, blunt boat hook-like tip  
**PAB-305.11**




**Kuglen Lens Manipulator and Iris Hook**  
Overall length 120 mm, ball-like tip  
**PAB-306.12**



**Kuglen Lens Manipulator and Iris Hook**  
Overall length 120 mm, angled ball-like tip  
**PAB-307.12**



**Kuglen Lens Manipulator and Iris Hook**  
Overall length 120 mm,  
**PAB-308.12** Straight  
**PAB-309.12** Angled




**Bonn Micro Iris Hook**

Overall length 120 mm very delicate 0,3 mm blunt hook

**PAB-310.12** Straight

**PAB-310.22** Angled




**Knolle Iris Hook**

Overall length 120 mm, delicate 0,3 mm blunt hook

**PAB-311.12** Straight

**PAB-312.12** Angled




**Shepard Iris Hook**

Overall length 120 mm small 1,0 mm blunt hook

**PAB-313.12** Straight


**PAB-313.22** Angled



**Hirschmann Iris Hook**

Overall length 120 mm, blunt 0,75 mm hook


**PAB-314.12**



**Osher Iris Hook**

Overall length 120 mm, blunt 1,5 mm hook

**PAB-316.12**




**Greather Multi-function Instrument**

Overall length 125 mm

**PAB-317.12** Straight


**PAB-317.22** Angled




**Osher Micro Suture Hook**

Overall length 120 mm


**PAB-318.12**




**Nordan Capsule Elevator**  
Overall length 115 mm, angled shaft with delicate ring  
**PAB-319.11**




**Anis Capsule Scrubber**  
Overall length 115 mm, with ball shaped tip  
**PAB-320.11**




**Anis Capsule Scrubber**  
Overall length 115 mm, with disc shaped tip  
**PAB-321.11**




**Mackool Capsule Retractor**  
Overall length 120 mm, 1,0 mm blunt hook with suspend 0,5 mm from end of tip  
**PAB-322.12**




**Maro Polishing Ring Curette**  
Overall length 105 mm, 45 dg angled shaft with 1,5 mm Semi - sharp open ring  
**PAB-323.10**



**Maro Polishing Ring Curette**  
Overall length 105 mm, 90 dg angled shaft with 1,5 mm semi-sharp open ring  
**PAB-324.10**



**Maro Polishing Ring Curette**  
Overall length 105 mm, 180 dg backward angled shaft with 1,5 mm semi-sharp open ring  
**PAB-326.10**




**Gills-Volker Capsule Polishing Curette**  
Overall length 105 mm, gently backward angled shaft with Ø1,5 mm sharp cup  
**PAB-327.10**



- Eye surgery • Augen-Instrumentarium • Oftalmologia
- Chirurgie ophtalmologique • Oftalmologia





**Straight Dressing Forceps with serration**

Overall length 100 mm



**PAB-330.05** 0,5 mm wide tips

**PAB-330.07** 0,7 mm wide tips

**PAB-330.10** 1,0 mm wide tips

**PAB-330.12** 1,2 mm wide tips

**PAB-330.14** 1,4 mm wide tips



**Straight Tissue Forceps with 1 x 2 teeth**

Overall length 100 mm


**PAB-332.05** 0,5 mm wide tips

**PAB-332.07** 0,7 mm wide tips

**PAB-332.10** 1,0 mm wide tips

**PAB-332.12** 1,2 mm wide tips


**PAB-332.14** 1,4 mm wide tips



**Bishop Harmon Dressing Forceps**

Overall length 85 mm criss-cross serrated

**PAB-333.06** 0,6 mm wide tips




**Bishop Harmon Tissue Forceps**

Overall length 85 mm 1 x 2 teeth

**PAB-334.06** Delicate 0,6 mm wide tips



**Gently curved Dressing Forceps with serration**  
 Overall length 100 mm gently curved tips  
**PAB-335.07** 0,7 mm wide tips




**Gently curved Tissue Forceps with 1 x 2 teeth**  
 Overall length 100 mm gently curved tips  
**PAB-336.07** 0,7 mm wide tips




**Strong Curved Dressing Forceps with serration**  
 Overall length 100 mm strong curved tips  
**PAB-337.07** 0,7 mm wide tips



**Strong curved Tissue Forceps with 1 x 2 teeth**  
 Overall length 100 mm Strong curved tips  
**PAB-338.07** 0,7 mm wide tips



**Angled Hess Iris Forceps Swiss pattern**  
 1 x 2 teeth, overall length 100 mm  
**PAB-339.10**




**Straight Graefe Iris Forceps**

Fine serrated, overall length 65 mm

**PAB-340.05** 0,5 mm wide shanks

**PAB-340.07** 0,7 mm wide shanks

**PAB-340.10** 1,0 mm wide shanks




**Straight Graefe Tissue Forceps**

1 x 2 teeth, overall length 65 mm

**PAB-341.05** 0,5 mm wide shanks

**PAB-341.07** 0,7 mm wide shanks

**PAB-341.10** 1,0 mm wide shanks




**Curved Graefe Iris Forceps**

1 x 2 teeth, overall length 65 mm

**PAB-342.05** 0,5 mm wide shanks

**PAB-342.07** 0,7 mm wide shanks

**PAB-342.10** 1,0 mm wide shanks




**Curved Graefe Tissue Forceps**

1 x 2 teeth, overall length 65 mm

**PAB-344.05** 0,5 mm wide shanks

**PAB-344.07** 0,7 mm wide shanks

**PAB-344.10** 1,0 mm wide shanks




**Straight Bonn Micro Iris Forceps**

0,12 mm 1x2 teeth, without tying platform

**PAB-345.75** Overall length 75 mm

**PAB-345.95** Overall length 95 mm





**Angled Bonn Micro Iris Forceps**

Angled 40 dg 0,12 mm 1 x 2 teeth, without tying platform

Overall length 72 mm

**PAB-346.72**




**Straight Bonn Micro Iris Forceps**

Ultra fine shaft with micro serrations,

Overall length

72 mm

**PAB-347.72**




**Angled Bonn Micro Iris Forceps**

Angled 40 dg ultra fine shaft with micro serrations,

Overall length 72 mm


**PAB-348.72**



**Strong Curved Graefe Tissue Forceps**

1 x 2 teeth, overall length 75 mm


**PAB-349.75**



**Straight Iris Tissue Forceps**

1 x 2 teeth, overall length 75 mm


**PAB-350.75**




**Curved Iris Forceps Swiss pattern**

1 x 2 teeth, overall length 100 mm


**PAB-351.10**




**Hess-Engelking Iris Forceps angled**  
 Overall length 72 mm  
**PAB-352.72** 1 x 2 teeth



**Gill Iris Forceps curved**  
 Overall length 75 mm  
**PAB-353.10** Cross - serrated  
**PAB-353.11** 1 x 2 teeth



**Iris Forceps Zuerich pattern**  
 Overall length 70 mm Extremely delicate  
**PAB-354.70** 1 x 2 teeth



**Botvin Iris Forceps double curved**  
 Overall length 70 mm  
**PAB-355.03** Serrated  
**PAB-355.04** 1 x 2 teeth



**Graefe Fixation Forceps without lock**

Overall length 110 mm

**PAB-356.35** Narrow, 3,5 mm jaws


**PAB-356.45** Wide, 4,5 mm jaws



**Graefe Fixation Forceps with lock**

Overall length 110 mm

**PAB-357.45** Wide 4,5 mm teeth



**Jess Fixation Forceps**

Overall length 110 mm

**PAB-358.11**



**Barraquer Fixation Forceps**

Overall length 100 mm


**PAB-359.10**




**Barraquer-Llobera Fixation Forceps**

Overall length 95 mm


**PAB-360.95**




**Green Fixation Forceps**  
 Overall length 100 mm Multiple teeth  
**PAB-362.05** 5,0 mm  
**PAB-362.10** 10,0 mm




**Elschnig Fixation Forceps**  
 Overall length 105 mm 1 x 2 angular teeth  
**PAB-363.03** 0,3 mm  
**PAB-363.09** 0,9 mm




**Lester Fixation Forceps straight**  
 Overall length 95 mm 0,6 mm teeth  
**PAB-364.12** 1 x 2 teeth  
**PAB-364.23** 2 x 3 teeth




**Kuhnt Fixation Forceps**  
 Overall length 110 mm straight shafts with curved  
 1 x 2 teeth  
**PAB-366.11**




**Bracken Scleral Fixation Forceps**  
 Overall length 105 mm straight shafts with  
 delicate 0,7 mm, 1x2 teeth  
**PAB-368.10**




**Farkas Bracken Corneal Fixation Forceps**  
 Overall length 105 mm straight shaft with delicate  
 0,12mm, 1x2 teeth  
**PAB-370.10**




**Troutman Superior Rectus Forceps**  
 Overall length 105 mm 45 dg angled shaft, 10 mm  
 from tip to bend, 1x2 teeth  
**PAB-372.10**




**Tennant-Troutman Superior Rectus Forceps**  
 Overall length 105 mm curved shanks 1x2 teeth  
**PAB-374.10**




**Sauer Suturing Forceps**  
 Overall length 90 mm straight shaft with 1x2  
 angled teeth  
**PAB-376.90**




**Manhattan Eye & Ear Forceps**  
 Overall length 90 mm straight shaft with 1x2 teeth  
 6mm long tying platform  
**PAB-378.90**




**Kremer Corneal Fixation Forceps straight**  
 Overall length 110 mm, with 2 point fixation, tips separated by 13mm 1x2 teeth each tip 0,12 mm, with catch lock  
**PAB-380.11**




**Kremer Corneal Fixation Forceps angled**  
 Overall length 110 mm, with 2 point fixation, tips separated by 13mm 1x2 teeth each tip 0,12mm, with catch lock  
**PAB-382.11**



**Bores Corneal Fixation Forceps Straight**  
 Overall length 103 mm, with 2 point fixation, tips separated by 3 mm 1x2 teeth each tip 0,12 mm  
**PAB-383.10**




**Rubman-Bores Corneal Fixation Forceps Curved.**  
 Overall length 103 mm, with 2 point fixation, tips separated by 3 mm 1x2 teeth each tip 0,12 mm  
**PAB-384.10**




**Arrowsmith Corneal Fixation Forceps Straight**

Overall length 107mm, with 2 point fixation, tips serrated by 8mm 1x2 teeth each side 0,12 mm with catch  
**PAB-386.10**




**Arrowsmith Corneal Fixation Forceps Curved**

Overall length 107mm, with 2 point fixation, tip serrated by 8mm 1x2 teeth each side 0,12 mm with catch  
**PAB-387.10**




**Bonaccolto Conjunctiva Forceps**

Overall length 110 mm with 16mm long platform serrated jaws and cross serrated tips  
**PAB-390.11** 1,2 mm wide  
**PAB-390.12** 1,7 mm wide  
**PAB-390.13** 2,4 mm wide




**McGregor Conjunctiva Forceps**

Overall length 110 mm with 5 mm long tying platform criss-cross serrated tips  
**PAB-392.11**



**Thorpe Conjunctiva Forceps**

Overall length 95 mm with tying platform and 2x3 teeth 0, set at 90 dg  
**PAB-394.95**




### Jaffe Tying Forceps

Overall length 100 mm, every delicate smooth jaws with 6mm long platform for 8-0 to 11-0 sutures

**PAB-396.10** Straight


**PAB-396.11** Curved



### Dardenne Tying Forceps angled

Overall length 100 mm extra delicate 7 mm long angled platform for 9-0 to 11-0 sutures

**PAB-398.10**




### Tennant Tying Forceps with Pin

Overall length 100 mm with guide pin and extra delicate 6 mm long platform for 9-0 to 11-0 sutures

**PAB-400.10** Straight

**PAB-400.11** Curved




### Harms-Tuebingen straight Tying Forceps

Overall length 100 mm with 6 mm long platform for 7-0 to 9-0 sutures

**PAB-402.01** Straight


**PAB-402.02** Curved



### Harms - Tuebingen Tying Forceps angled

Overall length 110mm with 7 mm long delicate angled platform for 7-0 to 9-0 sutures

**PAB-403.11**




**Troutman Tying Forceps**

Overall length 100 mm with extra fine 6 mm long platform  
for 7-0 to 9-0 sutures

**PAB-404.10** Straight


**PAB-404.11** Curved



**Bonn Tying Forceps (long model)**

Overall length 95mm with ultra fine 5 mm long platform  
for 9-0 to 11-0 sutures


**PAB-405.95**



**Bonn Tying Forceps (short model)**

Overall length 73 mm with ultra fine 5 mm long platform for  
9-0 to 11-0 sutures


**PAB-406.73**



**Micro Tying Forceps**

Overall length 85 mm with extra fine 6 mm long platform  
for 8-0 to 11-0 sutures


**PAB-407.85**



**Paufigue Tying Forceps**

Overall length 85 mm with standard 7 mm long platform sutures below 7-0


**PAB-408.85**



**Paufigue Tying Forceps**

Overall length 85 mm with delicate 5 mm long platform for 7-0 to 9-0 sutures


**PAB-410.85**



**Bechert-McPherson Tying and IOL Forceps**

Overall length 105 mm with angled 10 mm long platform for 8-0 to 11-0 sutures


**PAB-411.10**



**McPherson Tying Forceps angled**

Overall length 105 mm with 5 mm long platform for 7-0 to 9-0 sutures


**PAB-412.10**



**McPherson Tying Forceps Straight**

Overall length 105 mm with 5 mm long platform for 8-0 to 11-0 sutures


**PAB-413.10**



**Kelman-McPherson Tying Forceps**

Overall length 105 mm with angled 7,5 mm long platform for 8-0 to 11-0 sutures

**PAB-414.10**




**Anis Tying Forceps**

Overall length 100 mm with extra delicate 7 mm long platform for 9-0 to 11-0 sutures

**PAB-416.10** Straight


**PAB-416.11** Curved



**Sinsky Tying Forceps**



Overall length 100 mm with very delicate 7 mm long platform for 9-0 to 11-0 sutures

**PAB-418.10** Straight




**Gaskin Fragment Forceps angled**

Overall length 110mm with guide pin and extra fine  
11 mm long angled platform  
**PAB-420.11**





**Kelman-McPherson Tying Forceps angled**

Overall length 110mm with guide pin and extra fine  
8 mm long platform for 9-0 to 11-0 sutures  
**PAB-421.11**




**Tennant curved Tying Forceps**

Overall length 110mm with guide pin and extra fine  
6 mm long platform for 9-0 to 11-0 sutures  
**PAB-422.11**



**Tennant straight Tying Forceps**

Overall length 110mm with guide pin and extra fine  
6 mm long platform for 9-0 to 11-0 sutures  
**PAB-423.10**




**Castroviejo Suturing Forceps**

Overall length 105 mm with tying platform and 1x2 teeth

**PAB-424.12** 0,12 mm

**PAB-424.30** 0,30 mm

**PAB-424.50** 0,50 mm



**Castroviejo Suturing Forceps (standard handle)**

Overall length 100 mm with tying platform and 1x2 teeth


**PAB-425.12** 0,12 mm

**PAB-425.30** 0,30 mm

**PAB-425.50** 0,50 mm

**PAB-425.90** 0,90 mm

**PAB-425.95** 1,50 mm





**Castroviejo Suturing Forceps (wide handle)**

Overall length 100 mm with tying platform and 1x2 teeth

**PAB-426.12** 0,12 mm

**PAB-426.30** 0,30 mm


**PAB-426.50** 0,50 mm





**Sinder-Castroviejo Suturing Forceps**

Overall length 110 mm with guide pin, tying platform and 1x2 teeth 0,12 mm



**PAB-427.11**




**Jaffe Suturing Forceps**  
Overall length 110 mm 1x2 teeth 0,12 mm set at 90 dg and very delicate platform  
**PAB-428.11**




**McPherson Suturing Forceps**  
Overall length 105 mm with 1x2 teeth and straight 5 mm long platform  
**PAB-429.10**





**McPherson Suturing Forceps**  
Overall length 105 mm with 1x2 teeth and straight 5 mm long platform  
**PAB-430.10**



**Kelman-McPherson Suturing Forceps**  
Overall length 105 mm with 8 mm long angled platform and 1x2 teeth 0,12 mm set at 90 dg  
**PAB-431.10**




**Hoffer-McPherson Suturing Forceps**  
Overall length 105 mm with 11 mm long angled platform and 1x2 teeth, 0,12 mm set at 90 dg  
**PAB-432.10**



**Troutman Suturing Forceps**

Overall length 100 mm with extra fine tying platform and 1x2 teeth, 0,12 mm set at 90 dg


**PAB-434.10**



**Troutman Suturing Forceps**

Overall length 100 mm with extra delicate tying platform and 1x2 teeth 0,12 mm set at 90 dg


**PAB-435.10**



**Bonn Suturing Forceps (long model)**

Overall length 95 mm with ultra fine tying platform and 1x2 teeth 0,12 mm set at 90 dg


**PAB-436.95**



**Bonn Suturing Forceps (short model)**

Overall length 73 mm with ultra fine tying platform and 1x2 teeth 0,12 mm set at 90 dg

**PAB-437.73**




**Micro Tying Forceps**

Overall length 85 mm with extra fine platform and 1x2 teeth 0,12 mm set at 45 dg

**PAB-438.12** 12 mm long platform


**PAB-438.13** 3 mm long platform



**Paufigue Universal Suturing Forceps**

Overall length 88 mm 1x2 teeth of 0,70 mm set at 45 dg and tying platform


**PAB-439.12**



**Paufigue Universal Suturing Forceps delicate**

Overall length 85 mm 1x2 teeth of 0,30 mm set at 45 dg and tying platform

**PAB-440.85**




**Castroviejo Suturing Forceps (standard handle)**

Overall length 105 mm 1x2 teeth and angled tying platform

**PAB-441.30** 0,30 mm


**PAB-441.50** 0,50 mm

**PAB-441.55** 1,50 mm




**Castroviejo Suturing Forceps**

Overall length 110 mm with guide pin, tying platform and 1x2 teeth 0,12 mm set at 45 dg  
**PAB-442.11**




**Bonn Suturing Forceps**

Overall length 110 mm with guide pin, tying platform and 1x2 teeth 0,12 mm set at 90 dg  
**PAB-443.11**




**Kelman Suturing Forceps angled**

Overall length 110 mm with 8mm long tying platform, guide pin and 1x2 teeth, 0,12 mm set at 90 dg  
**PAB-445.11**




**Hoffer Suturing Forceps angled**

Overall length 110 mm with 11 mm long tying platform, guide pin and 1x2 teeth, 0,12 mm set at 90 dg  
**PAB-446.11**



**Troutman Suturing Forceps**

Overall length 110 mm with guide pin, tying platform and 1x2 teeth 0,12 mm set at 45 dg  
**PAB-447.11**




**Pierse Corneal Forceps**

Overall length 110 mm 1x1 pierse teeth 0,1 mm

**PAB-448.01** Straight

**PAB-448.02** Curved




**Pierse Corneal Fixation Forceps**

Overall length 75 mm 1x1 pierse teeth and tying platform

**PAB-448.10** 0,10 mm teeth

**PAB-448.20** 0,20 mm teeth



**Pierse - Polack Corneal Fixation Forceps**

Overall length 75 mm with 1x1 pierse teeth 0,20 mm at each tip, 1 mm spread


**PAB-449.75**



**Polack Corneal Fixation Forceps**

Overall length 75 mm, two 2,75 mm long tips separated 1 mm with 1x2 teeth 0, 12 mm


**PAB-450.75**



**Standard Corneal Suturing Forceps**

Overall length 72 mm with 1x2 teeth of 0,45 mm and tying platform

**PAB-451.72**





**Barraquer Corneal Utility Forceps**

Overall length 75 mm 1x2 teeth without tying platform

**PAB-452.12** 0,12 mm

**PAB-452.30** 0,30 mm




**Barraquer Corneal Suturing Forceps**

Overall length 75 mm 1x2 teeth with tying platform

**PAB-453.12** 0,12 mm


**PAB-453.30** 0,30 mm



**Barraquer Suture Removing Forceps**

Overall length 75 mm 1x1 teeth with tying platform

**PAB-454.75**




**Troutman Corneal Fixation Forceps**

Overall length 75 mm 1x2 teeth without tying platform

**PAB-455.12** 0,12 mm

**PAB-455.30** 0,30 mm




**Troutman Corneal Suturing Forceps**

Overall length 75 mm 1x2 teeth with tying platform

**PAB-456.12** 0,12 mm


**PAB-456.30** 0,30 mm



**Straight Corneal Tying Forceps**

Overall length 85 mm with extra fine 5 mm long jaws and tying platform



**PAB-457.85**



**Curved Corneal Tying Forceps**

Overall length 82 mm with extra fine 5 mm long jaws and tying platform


**PAB-458.82**



**Angled Corneal Tying Forceps**

Overall length 82 mm with extra fine 7 mm long angled jaws and tying platform

**PAB-459.82**





**Colibri Corneal Forceps**



Overall length 83 mm with extra fine 5 mm long jaws and tying platform

**PAB-460.11** Smooth tips, for tying

**PAB-460.13** Pointed tips, for suture removal




**Straight Corneal Fixation Forceps**  
Overall length 85 mm 1x1 pierse teeth with tying platform  
**PAB-461.85** 0,10 mm teeth





**Curved Corneal Fixation Forceps**  
Overall length 82 mm 1x1 pierse teeth with tying platform  
**PAB-462.82** 0,10 mm teeth



**Angled Corneal Fixation Forceps**  
Overall length 82 mm 1x1 pierse teeth with tying platform  
**PAB-463.82** 0,10 mm teeth




**Colibri Corneal Fixation Forceps**  
Overall length 83 mm 1x1 pierse teeth with tying platform  
**PAB-464.83** 0,10 mm teeth



**Straight Corneal Suturing Forceps**

Overall length 85 mm with 1x2 teeth set at 45 dg and tying platform


**PAB-465.85** 0,12 mm teeth



**Curved Corneal Suturing Forceps**

Overall length 82 mm with 1x2 teeth set at 45 dg and tying platform


**PAB-466.82** 0,12 mm teeth



**Angled Corneal Suturing Forceps**

Overall length 82 mm with 1x2 teeth set at 45 dg and tying platform



**PAB-467.82** 0,12 mm teeth



**Colibri Corneal Suturing Forceps**

Overall length 82 mm with 1x2 teeth set at 45 dg and tying platform


**PAB-468.82** 0,12 mm teeth



**Straight Corneal Utility Forceps**

Overall length 85 mm 1x2 teeth set at 90 dg with tying platform


**PAB-469.85** 0,12 mm teeth



**Curved Corneal Utility Forceps**

Overall length 82 mm 1x2 teeth set at 90 dg with tying platform


**PAB-470.82** 0,12 mm teeth



**Angled Corneal Utility Forceps**

Overall length 82 mm 1x2 teeth set at 90 dg with tying platform


**PAB-471.82** 0,12 mm teeth



**Colibri Corneal Utility Forceps**

Overall length 83 mm 1x2 teeth set at 90 dg with tying platform

**PAB-472.83** 0,12 mm teeth




**Girard Corneoscleral Forceps**

Overall length 100 mm 1x2 teeth 0,12 mm with tying platform

**PAB-473.01R** Right handed

**PAB-473.02L** Left handed



**Colibri Corneoscleral Forceps**

Overall length 100 mm with tying platform

**PAB-474.10** 1x2 teeth of 0,12 mm set at 45 dg

**PAB-474.12** 1x1 pierse teeth of 0,10 mm




**Pierse Corneoscleral Forceps**

Overall length 100 mm with tying platform and 1x1 pierse teeth 0,10 mm

**PAB-475.01** Straight

**PAB-475.02** Curved




**Suturing Corneoscleral Forceps**

Overall length 100 mm with 1x2 teeth of 0,12 mm and tying platform

**PAB-476.01** Straight

**PAB-476.02** Curved




**Tying Corneoscleral Forceps**


Overall length 100 mm with extra fine 5 mm long tying platform for 8-0 to 11-0 sutures

**PAB-477.01** Straight


**PAB-477.02** Curved




**Kuhnt Capsule Forceps**  
Overall length 100 mm 5x5 teeth  
**PAB-478.10**




**Swiss Model Capsule Forceps**  
Overall length 95 mm 6x6 teeth  
**PAB-479.95**




**Castroviejo Capsule Forceps**  
Overall length 100 mm, cross action  
**PAB-480.10**




**Arruga Capsule Forceps (new model)**  
Overall length 100 mm  
**PAB-481.10**



**Arruga Capsule Forceps**  
Overall length 100 mm  
**PAB-482.10**




**Capsulorhexis Forceps**  
 Overall length 100 mm extremely delicate 12 mm long vaulted shanks with ultra fine serrated tips  
**PAB-483.10**




**Lehner Capsulorhexis Forceps**  
 Overall length 105 mm with delicate 11 mm long vaulted shanks and ultra fine tips  
**PAB-484.10**





**Utrata Capsulorhexis Forceps**  
 Overall length 105 mm extra fine 11 mm long angled shanks with ultra fine tips  
**PAB-485.10**



**Nevyas Capsulorhexis Forceps**  
 Overall length 105 mm with extremely thin 11 mm long cystotome shaped tips for penetrating and grasping anterior capsule  
**PAB-486.10**




**Jaffe Capsulorhexis Forceps**  
 Overall length 105 mm with extremely thin 11 mm long shanks and extra delicate round grasping tips  
**PAB-487.10**




**Ernest - Calf Nucleus Cracking Forceps**

Overall length 105 mm extremely thin paddle shaped jaws with longitudinally serrated surfaces (cross action)  
**PAB-488.10**




**Kraff Nucleus Fracture Forceps**

Overall length 105 mm extremely thin paddle shaped jaws with cross serrated surfaces (cross action)  
**PAB-489.10**




**Kansas Nucleus Removal Forceps**

Overall length 105 mm with angled 10 mm long delicate jaws with multiple parallel teeth  
**PAB-490.10**



**Dardenne Nucleus Fragment Forceps**


Overall length 105 mm with angled 10 mm long shanks and 1 mm long grasping tips  
**PAB-491.10**



**Kelman-McPherson Fragment Forceps**

Overall length 100 mm with delicate angled 7.5 mm long shanks

**PAB-492.10**




**Gaskin Fragment Forceps**

Overall length 105 mm with delicate 45 dg angled shanks and smooth platform

**PAB-493.09** 9 mm angled


**PAB-493.11** 11 mm angled



**Osher Foreign Body Removal Forceps**

Overall length 105 mm extra fine 11 mm long 60 dg angled shanks with 3 mm long platform


**PAB-494.10**



**Kraff Capsule Tag Forceps**

Overall length 105 mm with thin 11 mm , long shanks and 1,5 mm long platform

**PAB-495.10**



**Witschel Capsule Flap Forceps**

Overall length 105 mm with 45 dg angled 11 mm long shanks and 5 mm long platform with 1x2 teeth of 0,10 mm

**PAB-496.10**



**Desmarres Chalazion Forceps**

Oval blades

**PAB-497.20** 20 mm

**PAB-497.26** 26 mm

**PAB-497.30** 30 mm



**Francis Chalazion Forceps**

Overall length 95 mm oval blades 12x14 mm


**PAB-498.95**



**Lambert Chalazion Forceps**

Overall length 90 mm round blade 8 mm dia

**PAB-499.90**



**Hunt Chalazion Forceps**

Overall length 95 mm round blade 12 mm dia

**PAB-500.95**



**Cauer Chalazion Forceps**

Overall length 100 mm round blade 12 mm serrated


**PAB-501.10**




**Lordan Chalazion Forceps**  
 Overall length 95 mm  
**PAB-502.95**




**Wies Chalazion Forceps**  
 Overall length 95 mm grooved upper jaw  
**PAB-503.95**



**Berke Ptosis Forceps**  
 Overall length 105 mm with longitudinally grooves and  
 slide lock  
**PAB-504.20** 20 mm  
**PAB-504.27** 27 mm



**Erhardt Lid Forceps**  
 Overall length 95 mm oval lower blade and serrated upper  
 blade  
**PAB-505.95**




**Snellen Entropium Forceps**  
 Overall length 95 mm  
**PAB-506.01** Right  
**PAB-506.02** Left



**Cueppers Strabismus Muscle Tucker**  
**PAB-507-01**



**Sattler Muscle Forceps**  
Overall length 100 mm, angled with 3 holes and teeth and side lock  
**PAB-508.22** Right  
**PAB-508.23** Left



**Jameson Muscle Forceps (child size)**  
Overall length 100 mm, angled with 4 holes and teeth and slide lock  
**PAB-509.01** Right  
**PAB-509.02** Left




**Jameson Muscle Forceps (adult size)**  
Overall length 100 mm, angled with 6 holes and teeth and slide lock  
**PAB-510.01** Right  
**PAB-510.02** Left



**Bangerter Muscle Forceps**  
Overall length 105 mm  
**PAB-511.10**


- Eye surgery • Augen-Instrumentarium • Oftalmologia
- Chirurgie ophtalmologique • Oftalmologia





**KALT**

**PAB-550.13**  
13cm  
gerade  
straight  
recta  
droit



**ARRUGA**

**PAB-551.13**  
13cm  
gebogen  
curved  
curva  
courbee

**PAB-552.14**  
14cm  
gerade  
straight  
recta  
droit


**PAB-553.14**  
14cm  
gebogen  
curved  
curva  
courbee



**CASTROVIEJO KALT**  
**PAB-555.16** 16CM




**STEVENS**  
**PAB-554.14** 14cm




Ohne Schloss  
without catch  
sin muelle sana  
ressort

**MICRO-CASTROVIEJO**  
PAB-560.09 9cm      PAB-561.09 9cm




Mit Schloss  
with catch  
con muelle  
avec ressort

**MICRO-CASTROVIEJO**  
PAB-562.09 9cm      PAB-563.09 9cm




Ohne Schloss  
without catch  
sin muelle  
sans resson

**CASTROVIEJO**  
PAB-564.14 14 cm      PAB-565.14 14 cm




Mit Schloss  
with catch  
con muelle  
avec ressort


**CASTROVIEJO**  
PAB-566.14 14 cm      PAB-567.14 14 cm



**MICRO - CASTROVIEJO**  
**PAB-570.14** 14cm  
**PAB-571.14** 14cm



**MICRO - CASTROVIEJO**  
**PAB-572.14** 14cm  
**PAB-573.14** 14cm



**CASTROVIEJO**  
**PAB-574.14** 14cm



**PAB-575.10**  
 10 cm

Ohne Schloss  
 without catch  
 sin muelle  
 sans ressort



**BARRAQUER**  
**PAB-576.13** 13cm  
**PAB-577.13** 13cm

Ohne Schloss  
 without catch  
 sin muelle  
 sans ressort



**BARRAQUER**  
**PAB-578.13** 13cm  
**PAB-579.13** 13cm

Mit Schioss  
 with catch  
 con muelle  
 avec ressort



Ohne Schloss  
without catch  
sin muelle  
sans ressort

**BARRAQUER**

**PAB-580.10**    **PAB-581.10**  
9.5 cm            9.5 cm



Mit Schloss  
with catch con  
muelle avec ressort

**BARRAQUER**

**PAB-582.10**    **PAB-583.10**  
9.5 cm            9.5 cm



Mit Schloss  
with catch  
con muelle  
avec ressort

**PAB-584.14**  
13.5 cm



Mit Schloss  
with catch  
con muelle  
avec ressort

**BARRAQUER**  
**PAB-585.13**  
13 cm

0.5mm



0.75mm



Ohne Schloss  
without catch  
sin muelle  
sans ressort

0.5mm



0.75mm



mit Schloss  
with catch  
con muelle  
avec ressort

**BARRAQUER**  
PAB-587.12 PAB-588.12  
11.5 cm

**BARRAQUER**  
PAB-589.12 PAB-590.12  
11.5 cm

0.75 mm



1.0 mm



Ohne Schloss  
without catch  
sin muelle  
sans ressort

0.75 mm



1.0 mm



mit Schloss  
with catch  
con muelle  
avec ressort


**BARRAQUER**  
PAB-591.14 PAB-592.14  
14 cm

**BARRAQUER**  
PAB-593.14 PAB-594.14  
14 cm




- Eye surgery
- Augen-Instrumentarium
- Oftalmologia
- Chirurgie ophtalmologique
- Oftalmologia






**PAB-595.11** 10,5 cm  
Gerade  
straight  
recta  
droit




**PAB-596.11**  
Gebogen  
curved  
curva  
courbe




**PAB-597.10** 9,5 cm  
Gerade  
straight  
recta  
droit


**PAB-598.10**  
Gebogen  
curved  
curva  
courbe




**STEVENS**  
**PAB-599.10** 10 cm




**STEVENS**  
**PAB-600.10** 10 cm




**STEVENS**  
**PAB-601.10** 10 cm




**STEVENS**  
**PAB-602.10** 10 cm




**BONN**  
**PAB-603.09**  
9 cm




**BONN**  
**PAB-604.09**  
9 cm




**MICROSCOPIC**  
**PAB-605.11**  
11 cm




**MICROSCOPIC**  
**PAB-606.11**  
11 cm




**PAB-607.10**  
9,5 cm




**PAB-608.10**  
9,5 cm



**STRABISMUS**  
**PAB-609.10**  
10 cm



**STRABISMUS**  
**PAB-610.10**  
10 cm




**PAB-612.09** 9.0 cm  
**PAB-612.11** 10.5cm  
**PAB-612.12** 11.5cm

**PAB-613.09** 9.0 cm  
**PAB-613.11** 10.5cm  
**PAB-613.12** 11.5cm

**PAB-614.12**  
 11.5 cm

**PAB-615.12**  
 11.5 cm

**PAB-616.12**  
 11.5 cm




**KNAPP**  
**PAB-617.11** 10.5cm  
**PAB-617.12** 12.0cm


**KNAPP**  
**PAB-618.11** 10.5cm  
**PAB-618.12** 12.0cm

**KNAPP**  
**PAB-619.11** 10.5cm  
**PAB-619.12** 12.0cm


**KNAPP**  
**PAB-620.11** 10.5cm  
**PAB-620.12** 12.0cm




**SREVENS**  
**PAB-622.11** 10,5 cm  
**PAB-622.12** 11,5 cm



**STEVENS**  
**PAB-623.11** 10,5 cm  
**PAB-623.12** 11,5 cm



**STEVENS**  
**PAB-624.11** 10,5 cm  
**PAB-624.12** 11,5 cm




**STEVENS**  
**PAB-625.11** 10,5 cm  
**PAB-625.12** 11,5 cm




**STRABISMUS**  
**PAB-626.11**  
 11 cm




**STRABISMUS**  
**PAB-627.11**  
 11 cm




**LAGRANGE**  
**PAB-628.09**  
 9 cm




**PAB-629.02**



**PAB-629.03**




**PAB-629.04**




**PAB-630.12** 11,5 cm  
**Fur Kinder**  
**for children**  
**para ninos**  
**pour enfants**


**PAB-629.01** 12 cm  
**PAB-629.02** 12 cm  
**PAB-629.03** 12 cm  
**PAB-629.04** 12 cm



**PAB-631.01**



**PAB-631.02**



**PAB-631.03**



**PAB-631.01** 13,5 cm  
**PAB-631.02** 13,5 cm  
**PAB-631.03** 13,5 cm

**WERBS DCR SCISSORS**  
 11 cm, 30 mm blades,  
 with 5mm angled ends



**PAB-635.02**



**PAB-636.02**



**PAB-637.02**



**PAB-638.02**



**PAB-639.02**



**PAB-640.02**



**PAB-635.03**



**PAB-636.03**



**PAB-637.03**



**PAB-638.03**



**PAB-639.03**



**PAB-640.03**



**PAB-635.04**



**PAB-636.04**



**PAB-637.04**



**PAB-638.04**



**PAB-639.04**



**PAB-640.04**



**PAB-635.02** 10cm  
**PAB-635.03** 10cm  
**PAB-635.04** 10cm

**PAB-636.02** 10cm  
**PAB-636.03** 10cm  
**PAB-636.04** 10cm

**PAB-637.02** 11cm  
**PAB-637.03** 11cm  
**PAB-637.04** 11cm

**PAB-638.02** 11cm  
**PAB-638.03** 11cm  
**PAB-638.04** 11cm

**NOYES**  
**PAB-639.02** 11.5cm  
**PAB-639.03** 11.5cm  
**PAB-639.04** 11.5cm

**NOYES**  
**PAB-640.02** 11.5cm  
**PAB-640.03** 11.5cm  
**PAB-640.04** 11.5cm



**TUEBINGEN**  
**PAB-642.85**  
8.5 cm



**TUEBINGEN**  
**PAB-643.85**  
8.5 cm




**PAB-644.80** 8.5 cm




**PAB-645.80** 8.5 cm

Trabekular-Schere  
Trabekular scissors  
Tijeras trabeculares  
Ciseaux a trabeculaires

Sharp pointed tips.  
11 mm from mid screw  
to tip



**CASTROVIEJO -  
VANNAS**  
**PAB-646.12**  
11.5 cm



**CASTROVIEJO -  
VANNAS**  
**PAB-647.12**  
11.5 cm



**McPHERSON-VANNAS**  
**PAB-648.85**  
8,5 cm



**VANNAS**  
PAB-649.80  
8 cm



**VANNAS**  
PAB-650.80  
8 cm



**VANNAS**  
PAB-651.80 8 cm  
PAB-652.80 8 cm



**VANNAS**  
PAB-652.80  
8 cm


Extra klein - extremely delicate - delicate - extra fine



**VANNAS**  
PAB-653.80  
8 cm



**VANNAS**  
PAB-654.80  
8 cm



**VANNAS**  
PAB-655.80  
8 cm



**VANNAS**  
PAB-656.80  
8 cm



PAB-662.01



PAB-662.02



PAB-663.01



PAB-663.02



PAB-662.03



PAB-662.04



PAB-663.03



PAB-663.04



PAB-658.55



PAB-660.55



PAB-659.55




PAB-661.55




**BARRAQUER**  
PAB-658.55 5.5 cm  
PAB-659.55 5.5 cm



**BARRAQUER**  
PAB-660.55 5.5 cm  
PAB-661.55 5.5 cm



**WECKER**  
PAB-662.01 11.5cm  
PAB-662.02 11.5cm  
PAB-662.03 11.5cm  
PAB-662.04 11.5cm




**WECKER**  
PAB-663.01 11.5cm  
PAB-663.02 11.5cm  
PAB-663.03 11.5cm  
PAB-663.04 11.5cm




**GIL**  
**PAB-664.01**




**ONG**  
**PAB-665.02**




**UTHOFF**  
**PAB-666.03**




**URIBE-STERN**  
**PAB-667.04**



**HOFFER-STERN**  
**PAB-668.05**




**PAB-670.10**  
10cm



**PAB-671.10**  
10cm



**PAB-672.10**  
10cm



**PAB-673.10**  
10cm



Kniegebogen  
angled to side  
anguloso  
coude

**PAB-674.10**  
10cm



**PAB-675.10**  
10cm



Aufgewinkelt  
angled on flat  
anguloso  
angulaire

**PAB-676.10**  
10cm




**PAB-677.10**  
10cm




Links  
left  
izquierda  
a gauche

**BARRAQUER  
PAB-678.01**




Rechts  
right  
derecha  
a droite

**BARRAQUER  
PAB-679.01**




Links  
left  
izquierda  
a gauche

**MC QUIRE  
PAB-680.01**




Rechts  
right  
derecha  
a droite

**MC QUIRE  
PAB-681.01**




Rechts  
right  
derecha  
a droite

**TROUTMAN-  
CASTROVIEJO  
PAB-682.01**




Links  
left  
izquierda  
a gauche

**TROUTMAN-  
CASTROVIEJO  
PAB-683.01**




Links  
left  
izquierda  
a gauche

**TROUTMAN  
PAB-684.01**




Rechts  
right  
derecha  
a droite

**TROUTMAN  
PAB-685.01**




Links  
left  
izquierda  
a gauche

**KATZIN  
PAB-686.01**




Rechts  
right  
derecha  
a droite

**KATZIN  
PAB-687.01**




Links  
left  
izquierda  
a gauche

**KATZIN  
PAB-688.01**



Rechts  
right  
derecha  
a droite

**KATZIN  
PAB-688.01**



Rechts  
right  
derecha  
a droite

**KATZIN  
PAB-689.01**



Links  
left  
izquierda  
a gauche

**KATZIN  
PAB-690.01**



Rechts  
right  
derecha  
a droite

**CASTROVIEJO  
PAB-691.10**  
10cm



Links  
left  
izquierda  
a gauche

**CASTROVIEJO  
PAB-692.10**  
10cm




Rechts  
right  
derecha  
a droite

**MINI-CASTROVIEJO  
PAB-693.10**  
10cm



Links  
left  
izquierda  
a gauche

**MINI-CASTROVIEJO  
PAB-694.10**  
10cm




Rechts  
right  
derecha  
a droite

**MINI-CASTROVIEJO  
PAB-695.10**  
10cm




Links  
left  
izquierda  
a gauche


**MINI-CASTROVIEJO  
PAB-696.10**  
10cm



**MINI-CASTROVIEJO  
PAB-697.11**  
10.5cm




**MINI-CASTROVIEJO  
PAB-698.11**  
10.5cm




**MICRO-CASTROVIEJO**  
**PAB-699.10**  
10cm



**MICRO-CASTROVIEJO**  
**PAB-700.10**  
10cm




**MINI-CASTROVIEJO**  
**PAB-702.13**  
12.5cm



**MINI-CASTROVIEJO**  
**PAB-703.13**  
12.5cm



**CASTROVIEJO**  
**PAB-701.11**  
10.5cm



**CASTROVIEJO**  
**PAB-704.90**  
9 cm



**CASTROVIEJO**  
**PAB-705.90**  
9 cm




**CASTROVIEJO**  
**PAB-706.90** 9 cm  
**PAB-707.90** 9 cm



**PAB-706.90**



**PAB-707.90**



**CASTROVIEJO**  
**PAB-708.01** 10.5cm  
**PAB-708.02** 10.5cm  
**PAB-708.03** 10.5cm  
**PAB-708.04** 10.5cm



**PAB-708.01**




**PAB-708.02**



**PAB-708.03**




**PAB-708.04**



**WESTCOTT**  
**PAB-709.11**  
 11cm




**WESTCOTT**  
**PAB-710.11** 11cm  
**PAB-710.13** 13cm





**WESTCOTT**  
**PAB-711.11**  
 11cm



- Eye surgery
- Augen-Instrumentarium
- Oftalmologia
- Chirurgie ophtalmologique
- Oftalmologia





**CASTROVIEJO-SCHAECHER**  
 PAB-721.70 7 cm




**PAB-722.11**




**WILDER**  
 PAB-723.11  
 11 cm, 4 1/2"



**PAB-724.11**



**SCHIOETZ**  
 PAB-725.10



**ALLEN-SCHIOETZ**  
 PAB-726.11



**CASTROVIEJO**  
 PAB-727.12  
 12CM, 4 3/4"



**ZIEGLER**  
 PAB-728.12 1 + 2  
 PAB-728.34 3 + 4  
 PAB-728.56 5 + 6  
 14cm



**BOWMAN**  
**PAB-730.04** Fig. 0000/000  
**PAB-730.03** Fig. 000/00  
**PAB-730.02** Fig. 00/0  
**PAB-730.11** Fig. 0/1  
**PAB-730.12** Fig. 1/2  
**PAB-730.23** Fig. 2/3  
**PAB-730.34** Fig. 3/4  
**PAB-730.45** Fig. 4/5  
**PAB-730.56** Fig. 5/6  
**PAB-730.67** Fig. 6/7  
**PAB-730.78.** Fig. 7/8

Neusilber  
 German Silver  
 Metal blanco  
 Ruolz

**BOWMAN**  
**PAB-731.04** Fig. 0000/000  
**PAB-731.03** Fig. 000/00  
**PAB-731.02** Fig. 00/0  
**PAB-731.01** Fig. 0/1  
**PAB-731.12** Fig. 1/2  
**PAB-731.23** Fig. 2/3  
**PAB-731.34** Fig. 3/4  
**PAB-731.45** Fig. 4/5  
**PAB-731.56** Fig. 5/6  
**PAB-731.67** Fig. 6/7  
**PAB-731.78** Fig. 7/8

Silber 800  
 Steriling silver  
 plata argent



4/0 3/0 = 0,6 / 7,0 mm  
 3/0 2/0 = 0,7 / 0,8 mm  
 2/0 0 = 0,8 / 0,9 mm  
 0 / 1 = 0,9 / 1,0 mm  
 1 / 2 = 1,0 / 1,2 mm  
 2 / 3 = 1,2 / 1,3 mm  
 3 / 4 = 1,3 / 1,4 mm  
 4 / 5 = 1,4 / 1,5 mm  
 5 / 6 = 1,5 / 1,6 mm  
 6 / 7 = 1,6 / 1,8 mm  
 7 / 8 = 1,8 / 2,0 mm



**BOWMAN**  
**PAB-732.04** Fig. 0000/000  
**PAB-732.03** Fig. 000/00  
**PAB-732.02** Fig. 00/0  
**PAB-732.01** Fig. 0/1  
**PAB-732.12** Fig. 1/2  
**PAB-732.23** Fig. 2/3  
**PAB-732.34** Fig. 3/4  
**PAB-732.45** Fig. 4/5  
**PAB-732.56** Fig. 5/6  
**PAB-732.67** Fig. 6/7  
**PAB-732.78** Fig. 7/8

Neusilber  
 German Silver  
 Metal blanco  
 Ruolz

**BOWMAN**  
**PAB-733.04** Fig. 0000/000  
**PAB-733.03** Fig. 000/00  
**PAB-733.02** Fig. 00/0  
**PAB-733.01** Fig. 0/1  
**PAB-733.12** Fig. 1/2  
**PAB-733.23** Fig. 2/3  
**PAB-733.34** Fig. 3/4  
**PAB-733.45** Fig. 4/5  
**PAB-733.56** Fig. 5/6  
**PAB-733.67** Fig. 6/7  
**PAB-733.78** Fig. 7/8

Silber 800  
 Steriling silver  
 plata argent





**West Lacrimal Chisel**  
**PAB-734.16**  
 Straight 16cm





**West Lacrimal Gouge**  
**PAB-735.16**  
 16cm



**Freer Periosteal Elevator**  
**PAB-736.18**  
 18cm D/E sharp/blunt 4,5 mm wide



**Joseph Periosteal Elevator**  
**PAB-737.16**  
 16 cm 4 mm wide



**Tenzel D/E Periosteal Elevator**  
**PAB-738.19**  
 195 mm shaft marked at 10mm, 20mm,  
 and 30mm from tip of dissector end

• Eye surgery • Augen-Instrumentarium • Oftalmologia



Luftkanule  
 Air injection cannula  
 Cannula para aire  
 Canule a air

**PAB-750.04** 27G., 0.4MM  
**PAB-750.06** 23G., 0.6MM



Luftkanule  
 Air injection cannula  
 Cannula para aire  
 Canule a air

**PAB-751.04** 27g., 0.4MM  
**PAB-751.06** 23g., 0.6MM



**PAB-752.04** 27G., 0.4MM  
**PAB-752.06** 23G., 0.6MM



Luftkanule  
 Air injection cannula  
 Cannula para aire  
 Canule a air

**PAB-753.27** 27g  
**PAB-753.30** 30g



Healon injektionskanule  
 Healon injection cannula  
 Canula para inyeccion de healon  
 Canule pour injection d' healon

**PAB-754.01**



Traenenwegkanule  
 Lacrimal cannula  
 Canula para canal lacrimal  
 Canule a voie lacrymal

**ANEL**  
**PAB-755.01**



Traenenwegkanule  
 Lacrimal cannula  
 Canula para canal lacrimal  
 Canule a voie lacrymal

**ANEL**  
**PAB-756.01**




Traenenwegkanule  
 Lacrimal cannula  
 Canula para canal lacrimal  
 Canule a voie lacrymal

**PAB-757.01**



**PAB-758.01**



**PAB-760.01**



**PAB-761.01**



**PAB-762.01**



**PAB-763.01**

Tränenwegkanule, rostfrei  
 Lacrimal cannula, SS  
 canule lacrimal, inox  
 Canule a voies lacrymales, inox

Tränenwegkanule mit 2 seitlichen Öffnungen  
 Lacrimal cannula with 2 lateral holes  
 Canula lacrimal con 2 orificios lateral  
 Canule a voies lacrymales avec 2 orifices lateral



**MC INTYRE**  
**PAB-764.40**



**MC INTYRE**  
**PAB-765.21**



**ATKINSON**  
**PAB-766.40**



**TROUTMAN**  
**PAB-767.01**

Retrobulbare Kanule, 40x0.6mm  
 Retrobulbar needle, 40mm x 23g  
 Canule retrobulbar 40x0.6mm  
 Canule retrobulbaire, 40x0.6mm

Alpha - Chymotrypsin Kanule  
 Alpha - Chymotrypsin cannula  
 Alpha - Chymotrypsin canula  
 Alpha - Chymotrypsin canule

Kanule zur Cyclodialyse  
 Cyclodialysis cannula  
 Canula para ciclodialisis  
 Canule pour cyclodialyse




**PAB-768.10**



**PAB-769.15**



**PAB-770.11**



**PAB-771.12**



**PAB-772.31**

Kanule zur Cyclodialyse  
 Cyclodialysis cannula  
 Canula para ciclodialisis  
 Canule pour cyclodialyse

Kanule zur Cyclodialyse  
 Cyclodialysis cannula  
 Canula para ciclodialisis  
 Canule pour cyclodialyse

Offnung rechts  
 Right side opening  
 Apertura a la derecha  
 Ouverture a droite

Offnung links  
 Left side opening  
 Apertura a la izquierda  
 Ouverture a gauche

Spulkanule fur radiale keratotomie  
 Radial Keratotomy irrigating cannula  
 Canula para irrigacion para queratotomia radial  
 Canule irrigatoire poue la keratotomie radiale



**SHEPARD**  
**PAB-773.01**

Spulkanule, seitliche offnung  
 irrigating cannula, lateral opening  
 Canula para irrigacion, orofocio lateral  
 Canule irrigatoire, orifice lateral



**WEST**  
**PAB-774.20** 20 g  
**PAB-774.23** 23 g

Sondenkanule  
 Probe cannula  
 Canula tipo sonda  
 Canule - sonde



**BANGERTER**  
**PAB-775.46** 40 X 0.6mm  
**PAB-775.48** 40 X 0.8mm  
**PAB-775.66** 60 x 0.6mm  
**PAB-775.68** 60 x 0.8mm

Augenbrause  
 Irrigating anterior chamber cannula  
 Canula para irrigacion de camara anterior  
 Canule irrigatoire multitrous



**CRANZ**  
**PAB-776.90**

Kortex - Extractor  
 Cortex Extractor  
 Extractor para corteza  
 Extracteur de cortex



**MC INTYRE**  
**PAB-777.10** **PAB-778.11**  
 Gerade Gewinkelt  
 straight angled  
 recto anguloso  
 droit coudee



**MC INTYRE**  
**PAB-779.23** 23 g  
**PAB-779.26** 26 g  
**PAB-779.30** 30 g



**MC INTYRE**  
**PAB-780.23** 23 g  
**PAB-780.26** 26 g  
**PAB-780.30** 30 g



**KNOLLE**  
**PAB-781.23** 23 g  
**PAB-781.30** 30 g



**SIMCO**  
**PAB-782.13** 13 mm  
**PAB-782.21** 21mm

Spulkanule mit Irishakchen  
 Irrigating iris hook  
 Cancho oros de irrigacion  
 irrigateur avec crochet pour l'iris



**MC INTYRE**  
**PAB-783.20** **PAB-783.21**  
 Gerade Gewinkelt  
 straight angled  
 recto anguloso  
 droit coudee

• Chirurgie ophtalmologique • Oftalmologia

Kapselpolierer, diamantbeschichtet  
 Capsule polisher, diamond dusted  
 Pulidor para capsula, con diamante  
 Polisher capsulaire, diamante



**KRATZ**  
**PAB-785.01**

Kapselpolierer, diamantbeschichtet  
 Posterior capsule polisher, diamond dusted  
 Pulidor para capsula, con diamante  
 Polisher capsulaire, diamante



**RUBMAN - KATZIN**  
**PAB-786.01**



**JENSEN**  
**PAB-787.21**  
**PAB-787.25**

Kapselpolierer, diamantbeschichtet  
 Posterior capsule polisher, diamond dusted  
 Pulidor para capsula, con diamante  
 polissoir capsulaire, diamante



**KNOLLE**  
**PAB-788.23** 23 g  
**PAB-788.30** 30 g



**KNOLLE**  
**PAB-789.23** 23 g  
**PAB-789.30** 30 g



Saugkeile  
 Absorbing points  
 Cunas absorbentes  
 coins absorbents

Packung 10 Stuck, steril  
 Package of 10 ea., Sterile  
 Paquete de 10 piezas, esteril  
 Boite a 10 pieces, sterile

**PAB-790.00**



Handriff  
 Infusion handle  
 Mango para infusion  
 Manche pour perfusion

**MC INTYRE**  
**PAB-791.60**



Handriff  
 Infusion handle  
 Mango para infusion  
 Manche pour perfusion

**KANSAS**  
**PAB-792.70**



**Holth Corneoscleral Punch**

Overall length 87mm

**PAB-795.10** 1,0 mm wide

**PAB-790.15** 1,5 mm wide

**PAB-797.20,** 2,0 mm wide

# SURGICAL *REGION*

Manufacture & Exporters

Surgical Region provides precision-crafted instruments for surgical, dental, beauty, and manicure/pedicure professionals. From orthodontic pliers and dental kits to hemostats, speculums, and nail nippers, our tools are built for performance, durability, quality, and affordability supporting excellence in every procedure.



For Products Inquiries & Purchases,  
Please Contact us:

✉ [info@surgicalregion.com](mailto:info@surgicalregion.com)

☎ +92 345 0089786